

Two ways of processing sound, two different listening experiences

Research shows that when presented with two distinctly different sound designs, a substantial proportion of listeners show clear and consistent preferences for one type over the other.¹

~40%

of people have a clear preference between two different sound designs¹

What this means for hearing aid fitting

85%

of hearing care professionals agree that no single hearing aid technology is best for everyone²

When it comes to hearing, **it's not one-size-fits-all**. Understanding the fundamental differences between time-domain and frequency-domain processing may help audiologists to:

Recognise that **differences in sound experience** can arise from processing architecture, not only fitting parameters



Understand why two well-fitted devices may be **experienced differently** by patients




Appreciate that **individual variability** in sound preference can influence overall satisfaction



A systematic focus on sound preference broadens fitting options and supports **more individualized fittings**

Filter banks are the foundation of signal processing


Digital hearing aids process sound through filter banks that split the signal into multiple frequency bands.



Distinct sound qualities arise from the signal processing approaches enabled by different filter bank designs.

Time-domain processing

Uses filters modelled on the logarithmic organization of the human cochlea with narrower filters at low frequencies, where sensitivity is high, and broader filters at higher frequencies.



Variable-width filters mimic the human ear




Enables low processing delay

Suited to listeners who prefer authentic, more natural sound clarity

Frequency-domain processing

Analyses sound using relatively narrow, uniformly spaced filters, which enables detailed parallel spectral analysis and efficient processing.



Extra processing power



High-resolution data processing

Often favored by listeners who prefer enhanced speech contrast and foreground focus

Reference: 1. Balling LW et al. Auditory sound preference as a measurable dimension in hearing-aid fitting: Evidence from comparisons of time-domain and frequency-domain processing, February 2026.

Reference: 2. Olesen, 2025. Internal study of 150 HCPs in the US and Germany.